COD BO13416 100 mL

COD AC16434 500mL

CONCENTRATED WASHING SOLUTION**CONCENTRATED WASHING SOLUTION****INTENDED USE**

Solution for the prevention of contamination between reagents in BioSystems automated analyzers. For *in vitro* professional use only in the clinical laboratory.

COMPOSITION

W1. Solution: Non-ionic surfactant 100 g/L, NaOH 1 mol/L, preservatives, pH >11.6.

DANGER: H314: Causes severe skin burns and eye damage. H412: Harmful to aquatic life with long lasting effects. P280: Wear protective gloves/protective clothing/eye protection/face protection. P303+P361+P353: IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower. P305+P351+P338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. P310: Immediately call a POISON CENTER or doctor/physician.

STORAGE AND STABILITY

Store at 15-30°C.

Components are stable once opened until the expiry date marked in the label if kept at the recommended storage temperature, well closed and if contamination is prevented during use.

Once diluted, the solution is stable for 3 months.

Indications of deterioration: Presence of particulate material, turbidity.

WARNINGS AND PRECAUTIONS

Exercise the normal precautions required for handling all laboratory reagents. Safety data sheet available for professional user on request. Disposal of all waste material should be in accordance with local guidelines. Any serious incident that might occur in relation to the device shall be reported to BioSystems S.A.

REAGENT PREPARATION

Dilute the concentrated with distilled water in a proportion 1/200.

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